

# Galileo

## RDT Platform

The Galileo platform makes rapid testing easy and convenient by simplifying the test procedure with interlocked user steps that mitigate common user errors.



The AtomoRapid 'Galileo' RDT platform improves test performance and usability with an integrated safety lancet, unique blood collection and delivery system in an easy to use all-in-one device.

The interlocked user steps make the handheld platform suitable for professional point of care and self-test applications. It can be customised for most lateral flow test strips that use capillary blood samples.

PLATFORM FEATURES	BENEFITS
<b>All-in-one test device</b>	<ul style="list-style-type: none"> <li>• Greater convenience and usability</li> <li>• Integration of features minimises the need for training</li> <li>• Supports clinical waste control by reducing the number of contaminated and potentially hazardous test components for disposal</li> </ul>
<b>Interlocked, sequenced user steps</b>	<ul style="list-style-type: none"> <li>• Supports compliance with test procedures</li> <li>• Ensures fewer user errors</li> </ul>
<b>Contact activated, auto-retracting safety lancet</b>	<ul style="list-style-type: none"> <li>• Safeguards against needlestick injuries and cross-contamination</li> <li>• Supports compliance to best practice in capillary sampling</li> </ul>
<b>Innovative blood collection unit</b>	<ul style="list-style-type: none"> <li>• Simplifies test procedure and ensures easy blood collection</li> <li>• Improves blood volume accuracy</li> </ul>
<b>Controlled delivery of blood to the test strip</b>	<ul style="list-style-type: none"> <li>• Eliminates blood splatter by transferring blood sample cleanly</li> <li>• Delivers blood to the correct point on the test strip</li> </ul>

### PLATFORM SPECIFICATIONS

<b>Type</b>	Single use, disposable lateral flow RDT platform
<b>Device dimensions</b>	108mm x 29mm x 21mm
<b>Safety lancet needle gauge</b>	21G
<b>Needle penetration depth</b>	2.4mm
<b>Blood volume* collected</b>	Designed to accommodate 5µL to 20µL
<b>Test strip dimensions</b>	4 (±0.2)mm x 73.5 (±0.5)mm

\*other volumes may be possible, further testing will be required



Lance Finger



Collect Blood



Deliver Blood



Add Buffer



Read Result

## Award Winning Design

Atomo rapid diagnostic tests have been recognised internationally for innovative design, performance and ease-of-use.



Discover how partnering with Atomo Diagnostics can provide your business with an efficient and low-risk pathway from product development to commercialisation of your rapid diagnostic test

### About Atomo Diagnostics

Atomo is a specialist medical device design and original equipment manufacturing (OEM) company with unique and proprietary technologies that improve rapid diagnostic testing by emphasising ease-of-use, accuracy and performance.

Atomo's all-in-one RDT device platforms put the end user first and make it easy to screen for various chronic conditions and test for a range of infectious diseases.

### Quality Assurance & Regulatory

Atomo's operations are ISO 13485 certified for the design, development and manufacture of in vitro diagnostic medical devices.



### Contacts

#### Australia - Head Office

Level 2  
701-703 Parramatta Road,  
Leichhardt, NSW 2040

Telephone: +61 (0)2 9099 4750

[atomodiagnostics.com/oem/](http://atomodiagnostics.com/oem/)  
[info@atomodiagnostics.com](mailto:info@atomodiagnostics.com)

**atomo diagnostics**  
SIMPLY BETTER DIAGNOSTICS